Phasix[®] Mesh / Phasix[®] ST Mesh

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Martin DP, Williams SF. "Medical Applications of Poly-4-Hydroxybutyrate: a Strong Flexible Absorbable Biomaterial."^{#17}

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Odermatt EK, Funk L, Bargon R, Martin DP, Rizk S, Williams SF. "Monomax® Suture: A New Long-Term Absorbable Monofilament Suture Made from Poly-4-Hydroxybutyrate."^{#18}

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Albertsmeier M et al. "Evaluation of the safety and efficacy of Monomax® suture material for abdominal wall closure after primary midline laparotomy — a controlled prospective multicenter trial: ISSAAC."^{#19}

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
INTRODUCTION

The most common cause of abdominal surgery is abdominal wall hernia defects [1]. Each year, there are approximately twenty million hernia repair surgeries performed worldwide, of which almost 700,000 surgeries are performed in America [2]. While the anatomy of this highly prevalent disease is well-known, recurrence of hernias and post-surgical complications such as chronic pain, inflammation, and infection have been reported in over half the patients treated [2]. If left untreated, herniated organs can become incarcerated and blood deprived which may result in infection, sepsis, and ultimately death [1].

Surgical meshes used to treat abdominal wall hernia defects act as an artificial wall - providing mechanical strength and acting as a buffer that helps prevent recurrence of the hernia. In 1958, Dr. Usher introduced the first artificial *Marlex* mesh and since then over 70 different types of mesh materials have been tested and used for repairing abdominal wall hernia defects [3]. Currently, the most common treatment option is the polypropylene (PP) mesh [4]. However, PP is permanent and has been associated with many post-surgical complications due to isolation of mesh implant and the host tissue inflammatory response resulting in chronic inflammation and pain. Many studies have shown that biodegradable mesh made from poly-4-hydroxybutyrate (P4HB) has comparable repair strength and successful transfer of load bearing from mesh to the repaired abdominal wall has been demonstrated [5]. P4HB is a natural polymer that is produced by microorganisms for the purpose of regulating energy metabolism [6]. The P4HB mesh is a fully biodegradable mesh made from polymer produced by *Escherichia coli* K12 bacteria via transgenic fermentation techniques [6]. P4HB fully degrades in vivo within approximately 78 weeks into carbon dioxide and water, which are broken down by the body very quickly via the Krebs cycle and beta-oxidation [7].

There are three phases to the natural healing process that occurs after surgery or injury. The first phase is the inflammatory phase with incursion of macrophage precursors that starts at day 1, peaks at day 3 and persists until the healing process is complete [1]. The next phase is the connective tissue phase, which is followed by the differentiation phase when tension resistance reaches its highest levels [1]. Macrophages undergo specific differentiation depending on micro environmental signals received during these phases and they can either be classically activated into M1 type or the alternatively activated to M2 type macrophages. The M1 macrophages are pro-inflammatory while the M2 macrophages reduce inflammation, are immune regulators and they promote tissue remodeling. Determining whether a mesh material activates M1 or M2 type macrophages can help explain post-surgical complications. Previous studies have shown that biodegradable meshes prompt more favorable M2 cell response as compared to polypropylene mesh [7].

OBJECTIVE

The purpose of this pre-clinical study was to evaluate the host foreign body immune response to P4HB mesh and

compare it to that of the widely used polypropylene mesh. The main focus was studying the types of macrophages activated and obtaining quantitative data.

SUCCESS CRITERIA

The success criteria for this study was set at an approximate increase of 15% in the M2 immune cell activation up to 200µm away from the site of implantation which would significantly prevent inflammation and promote healing.

METHOD

For the in vivo part of this study, approximately 1.5 centimeter square of mesh (either P4HB or PP) was implanted in the abdominal wall thickness of Sprague-Dawley rat model. My role in this study was to compare the immune system response to each type of mesh material and to determine which mesh material promotes the anti-inflammatory immune cell response with a special focus on M1 and M2 type macrophages. Explanting the four post-surgical samples of each mesh material along with surrounding tissue for the various designated time points (3, 7, 14, 21, or 35 days) took approximately 3 days. The next step was to prepare the tissue samples obtained for staining by starting with tissue fixation using immersion. The explanted tissues were immersed in 10% formalin solution for 6 hours at room temperature. Next, the tissues were dehydrated by immersing them in 90% ethanol twice for thirty minutes each. The dehydrated tissues were then embedded in paraffin and approximately 5µm sections were cut and floated in a water bath. Finally, slides pre-coated with gelatin, were used to mount the tissue sections.

The next step was staining the tissues using antibodies specific for markers of macrophage types. For instance, the CD-68+ marker was used to identify the non-induced or Pan-Mφ macrophages. The M1 macrophages were identified using the CD-86+ marker and the CD-206+ marker was used to identify the M2 macrophages. This step was very important as it would help determine the immune cell response for each mesh material. However, the task of optimizing staining posed challenges due to the reaction of primary and secondary antibodies causing higher background fluorescence. In order to reduce background fluorescence, the dilution of proteins present in blocking solutions was key. This step helped minimize non-specific binding of antibodies.

The final step in this process was analyzing the stained tissue samples using a microscope. All the images taken were then uploaded to the Cell Profiler software so that the obtained data could be quantified. Cell quantification around the mesh material was done at four distances away from the mesh: 0-50 µm, 50-100 µm, 100-150 µm and 150-200 µm. After the data was quantified, we were able to compare the data for P4HB mesh to the control (Polypropylene mesh).

RESULTS

Significant differences were observed between the two mesh types. The absorbable P4HB mesh demonstrated a higher ratio of M2 over M1 macrophages especially with the earlier time points; day 3 and day 7, compared to polypropylene mesh.

This increase in the M2 (anti-inflammatory) response at the earlier time points is likely due to the increase in the expression of M2 related proteins such as *Ym1*, *Arg1*, and *Fizz1* and the secretion of cytokines that were polarized to M2 type in animals treated with P4HB mesh. Figure 1 shows how the M1 and M2 macrophages were quantified at the various distances from the site of implant. While the figure here shows the different types of immune cells around the implant, the main goal of this study was to quantify only the M1 and M2 type of macrophages.

Immune Cell Activation around Implant


Figure 1. Immune cell activation around the site of mesh implant at various distances. Ring 1 is 50 µm, Ring 2 is 100 µm, Ring 3 is 150 µm and Ring 4 is 200 µm away from mesh.

Figure 2 below shows the M2/M1 ratio for both mesh materials for different time points (Days 3, 7, 14, 21 and 35) at each of these distances (rings) from the mesh. The four bars on each graph for P4HB and PP represent data collected from explanted tissues of the four animals that were treated with each of these mesh types.


Figure 2. M1 vs. M2 macrophage activation at the site of implantation at Days 3,7,14,21 and 35. Standard error bars are presented. P4HB shows a higher M2 activation than PP.

The X-axis of the graphs show the mesh type that was tested while the Y axis represents the M2 vs M1 ratio of macrophages activated at that time point for each ring. The figure also shows standard error bars. The error bars represent the differences between these biological replicates. Preliminary results show a 33% increase in M2 type macrophages for P4HB mesh over Polypropylene mesh.

DISCUSSION

The results of this study showed that P4HB does stimulate more anti-inflammatory immune cell response as compared to polypropylene mesh. We had originally expected that there would be an approximate 15% increase in the M2 macrophage cell activation but the P4HB mesh performance far exceeded this expectation. As mentioned earlier, M2 macrophages are anti-inflammatory as opposed to M1 macrophages, which are pro-inflammatory. The results of this study echo the results of other previous studies that biodegradable meshes do in fact activate the M2 anti-inflammatory cells [4,5,7]. This can prove to be very useful in treating hernias as most post-surgical complications with hernia repairs are associated with inflammation at the site of implant. While polypropylene mesh is widely used for hernia repairs, the fact that almost 50% of the patients treated develop some post-surgical complications has highlighted the need for a better treatment option.

However, there are certain limitations to this study. The statistical analyses for cell quantification that need to be done are very complex. Other biodegradable mesh materials also need to be tested to compare performance. The anti-microbial qualities of the mesh materials should also be tested so that infections and sepsis can be prevented. This study focused on M1/M2 type immune cells, but other immune cells also need to be quantified and the results obtained in this study need to be confirmed with in vitro models as well. Despite some challenges and requirements for further testing, the results of many studies have shown that biodegradable meshes such as P4HB can in fact be the best potential solution for treating hernia defects.

ACKNOWLEDGMENTS

This work was supported as part of an Industrial Partnership with C.R. Bard Inc. I would also like to thank Denver Faulk, Department of Bioengineering, University of Pittsburgh and the Badylak Laboratory.

REFERENCES

- [1] Zogbi, L. (2011). *The Use of Biomaterials to Treat Abdominal Hernias* Biomaterials Applications for Nanomedicine <http://cdn.intechopen.com/pdfs-wm/23633.pdf>
- [2] Dabbas, N., Adams, K., Pearson, K., & Royle, G. (2011). *Frequency of abdominal wall hernias: is classical teaching out of date?* Journal of the Royal Society of Medicine; Short Reports, 2(1), 5. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3031184/>
- [3] Inguinal Hernia: Anatomy and Management http://www.medscape.org/viewarticle/420354_3
- [4] Pineda C., Sicari B.M., Londono R., Dziki J., Faulk D., Giglio R., Wang D., Gandhi R.M., Shaffiey S., Badylak S.F. (2015) *The Effect of Poly-4-Hydroxy Butyrate upon Host Innate Immune Response and Antimicrobial Activity* TERMIS poster (Manuscript in preparation)
- [5] Martin, D.P., Badhwar, A., Shah, D.V., Rizk, J.R., Scott, J.R. (2013) *Characterization of poly-4-hydroxybutyrate mesh for hernia repair applications* J Surg Res 184(2) <http://www.sciencedirect.com/science/article/pii/S0022480413002382>
- [6] Deeken, C. & Matthews, B. (2013) *Characterization of the Mechanical Strength, Resorption Properties, and Histologic Characteristics of a Fully Absorbable Material (P4HB-PHASIX Mesh) in a Porcine Model of Hernia Repair* ISRN Surg <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3679684/>
- [7] Martin D.P., Williams S.F. (2003) *Medical applications of poly-4-hydroxybutyrate: a strong flexible absorbable biomaterial.* *Biochemical Engineering Journal*;16 (2):97-105

Modulatory Mechanisms of Immune Host Response and Antimicrobial Peptide Expression induced by Poly(4-Hydroxy Butyrate)

Catalina Pineda^{1,2}, George S. Hussey^{1,3}, Ross M. Giglio¹, Angela Ramirez¹, Brian M. Sicari^{1,3}, Ricardo Londono⁴, Riddhi M. Gandhi^{1,2}, and Stephen F. Badylak, D.V.M., Ph.D., M.D.^{1,2,3}

¹McGowan Institute for Regenerative Medicine, ²Department of Bioengineering, ³Department of Surgery, ⁴School of Medicine, University of Pittsburgh, Pittsburgh, Pennsylvania.

Statement of Purpose: The bioactive metabolite, Butyrate, is a naturally occurring histone deacetylase (HDAC) inhibitor found mainly in the gastrointestinal tract. Butyrate has been associated with seemingly paradoxical anti-inflammatory and antimicrobial effects (Donohoe, 2012). In particular, butyrate promotes an anti-inflammatory M2-like macrophage activation (Ji, 2016) and expression of the antimicrobial peptides (AMPs) cathelicidins and β -defensins (Kida, 2006). Poly(4-hydroxy butyrate) (P4HB), a hydroxylated form of butyrate, is a byproduct of bacterial metabolism, that has been commercially produced and used as an FDA approved 'biosynthetic' material for multiple soft tissue surgical applications (Williams, 2013). Surgical mesh materials composed of P4HB degrade by hydrolysis *in vivo*, releasing the monomer γ -hydroxy butyrate (GHB), and have been shown to promote resistance to deliberate surgical site bacterial contamination (Pineda, 2016). The mechanism(s) behind the bacterial resistance remains only partially understood and is known to involve the host innate immune system. The present study evaluates the mechanisms by which a biosynthetic mesh material composed of P4HB modulates immune host response and promotes AMPs expression.

Methods: An *in vitro* accelerated hydrolysis method with hydrochloric acid (Martin, 2013) was used to obtain the degradation byproducts of P4HB. After further neutralization to physiological pH and dialysis to remove the excess of salt, the hydrolysis byproducts were used in mice bone marrow-derived macrophages to evaluate the activation of markers associated with pro-inflammatory (iNOS) and anti-inflammatory (Fizz-1) phenotype. The expression of the cathelicidin LL-37 and β -defensin-4 were also determined. Known factors promoters of pro-inflammatory (LPS and IFN- γ) or anti-inflammatory (IL-4) phenotypes, or commercial GHB were used as controls. In addition, surgical mesh materials composed of P4HB were evaluated *in vivo* using a rat bilateral partial thickness abdominal wall defect model (Sicari, 2013). Intermediaries of the metabolic pathway of the degradation products of P4HB were assessed by GC-MS. Immunolabeling quantification was used to determine the immunomodulatory effects of P4HB upon macrophage phenotype, using CD206 (anti-inflammatory), and CD-86 (pro-inflammatory) markers, and cathelicidin LL-37 expression, at 3, 7, 14, 21, and 35 days post-implantation. **Results:** The results suggest that surgical mesh materials composed of P4HB promote an anti-inflammatory M2 phenotype (Figures 1A & 2A) and increased expression of AMPs (Figures 1B-C & 2C). The metabolism of the degradation byproducts of P4HB increased the production to succinic acid (Figure 2B), an intermediary in the degradation of P4HB-derived monomers. Furthermore,

the mechanism by which P4HB influences the macrophages involves transcriptional regulation of AMPs through MAP kinases pathways and HDAC inhibition.


Figure 1. Modulatory effects of P4HB *in vitro*. **A.** Immunolabeling for F4/80, iNOS, and Fizz-1 markers. **B.** Quantitative RT-PCR of cathelicidin LL-37 and β -defensin-4, fold change relative to macrophages in MCSF control media. **C.** Immunolabeling of cathelicidin LL-37 and percentage of macrophages expressing the AMP.


Figure 2. Modulatory effects of P4HB *in vivo*. **A.** Macrophage phenotype surrounding the mesh implant. **B.** Production of succinic acid. **C.** Expression of cathelicidin LL-37 *in vivo*.

Conclusions: Bioscaffolds composed of P4HB have the potential to provide a constructive modulatory effect commonly seen with alternatively activated macrophages and have the advantage of a reproducibly manufactured "synthetic" bioscaffold.

References: (Donohoe DR, *et al.* Mol Cell. 2012;48:612-626) (Campbell Y, *et al.* Eur J Nutr. 2012;51:899-907) (Kida Y, *et al.* Mol Immunol. 2006;43:1972-1981) (Ji J, *et al.* Sci Rep. 2016;6:24838) (Williams SF, *et al.* Biomed Tech. 2013;58:439-452) (Pineda C, *et al.* Manuscript in Preparation. 2016) (Martin DP, *et al.* J Sur Res. 2013;184:766-773) (Sicari B, *et al.* Methods Mol Biol. 2013;1037:3-25).

Acknowledgements: This work was supported as part of an industrial partnership with CR Bard, Inc.

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


Prospective Evaluation of Poly-4-hydroxybutyrate Mesh in CDC Class I, High Risk Ventral and Incisional Hernia Repair: 3-Year Follow-up

Roth JS¹, Anthone GJ², Selzer DJ³, Poulouse BK⁴, Bittner JG⁵, Hope WW⁶, Dunn RM⁷, Martindale RG⁸, Goldblatt MI⁹, Earle DB¹⁰, Romanelli JR¹⁰, Mancini GJ¹¹, Greenberg JA¹², Linn JG¹³, Parra-Davila E¹⁴, Sandler BJ¹⁵, Deeken CR¹⁶, Voeller GR¹⁷

¹University of Kentucky (Lexington, KY), ²Methodist Health System (Omaha, NE), ³Indiana University (Indianapolis, IN), ⁴Vanderbilt University (Nashville, TN), ⁵Virginia Commonwealth University (Richmond, VA), ⁶New Hanover Regional Medical Center (Wilmington, NC), ⁷University of Massachusetts (Worcester, MA), ⁸Oregon Health & Science University (Portland, OR), ⁹The Medical College of Wisconsin (Milwaukee, WI), ¹⁰Baystate Medical Center (Springfield, MA), ¹¹University of Tennessee (Knoxville, TN), ¹²University of Wisconsin (Madison, WI), ¹³North Shore University (Evanston, IL), ¹⁴Celebration Health (Celebration, FL), ¹⁵University of California (San Diego, CA), ¹⁶Covalent Bio, LLC (St. Louis, MO), ¹⁷University of Tennessee (Memphis, TN)

Objective

Long-term resorbable mesh represents a promising technology for complex ventral and incisional hernia repair (VIHR). Preclinical studies indicate that poly-4-hydroxybutyrate (P4HB) bioresorbable mesh supports strength restoration of the abdominal wall. The objective of this prospective, multi-institutional study was to evaluate outcomes of CDC Class I/High-Risk patients undergoing VIHR with P4HB mesh.

Study Device

- Knitted monofilament P4HB (Phasix™ Mesh)**
- Predictably and gradually degrades via hydrolysis
- Contributes mechanical strength up to 12 months**
- Scaffold enables remodeling by host tissue**
- Essentially fully resorbed by 18 months**.


Figure 1: Phasix® Mesh

Inclusion Criteria

- CDC Class I wounds (clean)
- Hernia defect: 10-350 cm²
- ≤ 3 prior repairs
- ≥ 1 high-risk criteria (detailed in Table 1)

Table 1 — Incidence of high-risk comorbid conditions (%)

Comorbid Conditions	Percentage of Subjects (%)
BMI (30-40 kg/m ²)	78.5
Hypertension	59.5
Cardiovascular disease	34.7
Diabetes	33.1
COPD	28.1
Malignancy	24.8
Active smoker	23.1
Immunosuppression	8.3
Chronic corticosteroid use	5.0
Advanced age	5.0
Hypo-albuminemia	2.5
Renal insufficiency	0.8

Study Population

Table 2 — Subject demographics

Subjects treated	n=121
Subjects with 3-year-follow-up	n=82 (67.8 %)
Sex	n=46 (38 %) male; n=75 (62 %) female
Age (years)	54.7 ± 12.0
Body mass index (kg/m ²)	32.2 ± 4.5

Surgical Technique

Subjects underwent open ventral hernia repair with retrorectus or onlay placement of mesh (with or without additional component separation, CST). Per protocol, Phasix™ Mesh was to be positioned with the edges extending beyond the margins of the defect by least 5 cm. However, this was not achieved in several instances. Phasix™ Mesh was fixated either according to surgeon preferred technique or not fixated. In most cases, the hernia defect was closed by approximating the fascial edges, including additional CST, if required.

Table 3 — Operative details (CST: component separation)

Defect (cm ²), mean ± SD	115.7 ± 80.6
Mesh (cm ²), mean ± SD	459.38 ± 172.3
Surgical procedure time (hrs), mean ± SD	2.8 ± 1.4
Age (years)	
Retrorectus without CST, n (%)	43 (35.5 %)
Retrorectus with CST, n (%)	45 (37.2 %)
Onlay without CST, n (%)	24 (19.8 %)
Onlay with CST, n (%)	8 (6.6 %)
Other, n (%)	1 (0.8 %)

Data Collection

Postoperative physical exam, quality of life, pain VAS, adverse events, and pain medication usage were assessed at 1, 3, 6, 12, 18 and 36 months, with long-term (60 months) follow-up ongoing.

Long-Term Outcomes (3 Years)

Table 4 — Primary and representative secondary outcomes

Hernia recurrence	15.7 % (n=19)
Onlay	28.1 %*
Sublay	9.3 %*
Surgical site infection	9.1 % (n=11)
Seroma requiring intervention	6.6 % (n=8)
Rate of reoperation	11.6 % (n=14)
Device-related adverse events	15.7 % (n=19)

Conclusions

Intermediate term outcomes of VIHR with P4HB mesh are favorable and demonstrate low recurrence rates at the 3-year (36-month) postoperative time frame. 5-year (60-month) follow-up is ongoing.

Disclosures

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.

This study was sponsored by Davol Inc., a subsidiary of C. R. Bard, Inc. (Warwick, RI.) Bard and Davol have joined BD (Franklin Lakes, NJ).

Authors were reimbursed for expenses related to the conduct of the study.

JSR, GJA, JGB, WWH, RGM, MIG, DBE, GJM, JAG, EPD, BJS, CRD, and GRV are paid consultants for C. R. Bard, Inc. (Daval).

Presented at SAGES 2019 (Baltimore, MD)

* Presented at International Hernia Symposium 2019 (Jeekel H, Erasmus University Medical Center, Rotterdam)

** Preclinical data on file. Results may not correlate to clinical performance in humans.

A Post-Market, Prospective, Multi-Center, Single-Arm Clinical Investigation of Phasix Mesh for VHWG Grade 3 Midline Incisional Hernia Repair: A Study of Patients in the European Union

T. Tollens, MD¹, L.N. Jørgensen, MD, PhD², T.S. de Vries Reinlingh, MD, PhD³, F. Köckerling, MD, PhD⁴, G. Piessen, MD⁵, M. A.C.J. Windsor, MD, PhD⁶, M. Miserez, MD, PhD⁷, F. Berrevoet, MD, PhD⁸, D. Dousset, MD, PhD⁹, F. Gossetti, MD¹⁰, H.L. van Westreenen, MD¹¹, R. H. Fortelny, MD¹², G. Woeste, MD¹³, U. Pession, MD¹⁴, J.F. Lange, MD, PhD¹⁵, G. Tetteroo, MD¹⁵, A. Koch, MD¹⁶, J. Jeekel, MD, PhD¹⁵

¹Imelda Hospital (Bonheiden, Belgium); ²University of Copenhagen (Copenhagen, Denmark); ³Elkerliek Hospital (Helmond, The Netherlands); ⁴Vivantes Klinikum Spandau (Berlin, Germany); ⁵University Hospital Lille (Lille, France); ⁶University College London Hospital (London, United Kingdom); ⁷University Hospital Leuven (Leuven, Belgium); ⁸University Hospital Ghent (Ghent, Belgium); ⁹Hôpital Cochin (Paris, France); ¹⁰Università di Roma Sapienza (Rome, Italy); ¹¹Isala Zwolle (Zwolle, The Netherlands); ¹²Wilhelminenhospital (Vienna, Austria); ¹³Klinikum der Johann Wolfgang Goethe-Universität (Frankfurt am Main, Germany); ¹⁴University Hospital (Frankfurt, Germany); ¹⁵Erasmus University Medical Centre Rotterdam (Rotterdam, The Netherlands); ¹⁶Chirurgische Praxis Cottbus (Cottbus Area, Germany)

Background

Incisional hernia is a frequent complication of midline laparotomy. The use of mesh has been reported to lead to fewer recurrences compared to primary repair. However, in Ventral Hernia Working Group (VHWG) Grade 3 patients, whose hernia is potentially contaminated, synthetic mesh may be prone to infection.* There is a strong preference for biological mesh in contaminated fields, which may better resist infection with little chance of foreign body reaction compared to permanent synthetic mesh. However, non-crosslinked biological mesh products often degrade too quickly to facilitate native cellular ingrowth. Phasix[®] Mesh is a biosynthetic mesh with both the biocompatibility and resorbability of a biological mesh and the mechanical strength of a synthetic mesh.*

Objective

This multicenter, single-arm study aims to assess the use of Phasix[®] Mesh in VHWG Grade 3* European hernia patients.

Study Design

A total of n=96 VHWG Grade 3 hernia patients were enrolled in the trial, with n=84 treated with PPhasix[®] Mesh at 15 sites across Europe (Clinical Trials #NCT02720042). The primary outcome is Surgical Site Occurrence (SSO) requiring intervention, including hematoma, seroma, surgical site infection (SSI), wound dehiscence, skin necrosis, and fistula formation through 3 months. Secondary outcomes include hernia recurrence, infection, pain, and quality of life through 24 months. Follow-up visits are scheduled at drain removal, at discharge or staple removal, and at 1, 3, 6, 12, 18, and 24 months postoperatively.

Subjects enrolled in trial, n	96
Subjects implanted with Phasix [®] Mesh, n	84
Subjects with 12-month follow-up, n (%)	84 (100%)
Gender	
Male, n (%)	51 (60.7%)
Female, n (%)	33 (39.3%)
Race	
Asian, n (%)	1 (1.2%)
White / Caucasian, n (%)	72 (85.7%)
Not reported (France), n (%)	11 (13.1%)
Age (years), mean	62.5
Body mass index (kg / m ²), mean	27.8


Figure 1: Preoperative diagnoses where “multiply recurrent” encompasses second (6.0%), third (1.2%), and fourth (1.2%) time recurrences.

Results

As shown in Table 1, n=84 subjects have completed the 12-month follow-up visit. The majority are Caucasian (85.7%) males (60.7%) with a mean age of 62.5 years and mean BMI of 27.8 kg/m². Over 80% of the hernias were primary incisional hernias (Figure 1) in a variety of locations throughout the abdomen (Figure 2).

The hernia defect size was 109.2 ± 87.9 cm², and the surgical approach was primarily retrorectus (Figure 3) with or without Component Separation Technique (57.1% and 41.7%, respectively). SSOs requiring intervention at 12 months postimplantation are depicted in Figure 4, notably 2.4% skin necrosis, hematoma, and fistula (2 subjects each), 7.1% hernia recurrence (6 subjects), 8.3% seroma (7 subjects), 10.7% skin dehiscence (9 subjects), and 13.1% surgical site infection (11 subjects). One subject was Grade 4.

Conclusion

Phasix[®] Mesh demonstrated favorable outcomes in VHWG Grade 3 European patients, exhibiting reasonable incidence of SSO with low rate of hernia recurrences at 12 months.* Longer-term 24-month follow-up is ongoing.

Disclosures

This project was sponsored by Davol Inc., a subsidiary of C. R. Bard Inc. Bard and Davol have joined BD. (Franklin Lakes, NJ).

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions and instructions for use.

* No mesh is indicated for use in contaminated or infected fields.


Figure 2: Location of hernia where “Other” represents midline, umbilical + epigastric, subxiphoidal + umbilical, M2 + M3, supra-umbilical, multiple midline defects, umbilical + infraumbilical, midline + parastomal, and other combinations.


Figure 3: Surgical approach (CST = Component Separation Technique)


Figure 4: Primary endpoint — Surgical Site Occurrences (SSOs) requiring intervention at 3 months postimplantation.

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Rx Only

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Quellen

- Prospective Evaluation of Phasix[®] Mesh in CDC Class I/High Risk Ventral and Incisional Hernia Repair: 18 Months Follow-up, (Roth et al., 2017, Surgical Endoscopy)
- Complex and infected ventral hernia repair using P4HB meshes (Auer et al., 2017, Sages Kongress)
- Prospective Multicenter Evaluation of Ventral/Incisional Hernia Repair with P4HB (Novitsky et al., 2016, Sages Kongress)
- Initial Experience With Biologic Polymer Scaffold (Poly-4-hydroxybutyrate) in Complex Abdominal Wall Reconstruction (Buell et al., 2016, Annals of Surgery)
- A Post-Market, Prospective, Multi-Center, Single-Arm Clinical Investigation of Phasix[®] Mesh for VHWG Grade 3 Midline Hernia Repair (Jeekel et al., 2018, AHS/EHS Kongress Miami)
- Deeken CR, Abdo MS, Frisella MM, Matthews BD. "Physicomechanical evaluation of absorbable and nonabsorbable barrier composite meshes for laparoscopic ventral hernia repair." *Surgical Endoscopy* 25.5 (2010): 1541-552.
- Estimated from Standard Curve in (Martin DP, et al. "Characterization of poly-4-hydroxybutyrate mesh for hernia repair applications." *Journal of Surgical Research*. 2013; (184): 766-773.
- Characterization of host response, resorption, and strength properties, and performance in the presence of bacteria for fully absorbable biomaterials for soft tissue repair, Stoikes et al., 2017
- Evaluation of the safety and efficacy of Monomax[®] suture material for abdominal wall closure after primary midline laparotomy — a controlled prospective multi-centre trial: ISSAAC, Arbeitsmeier et al., 2011
- Phasix[®] Mesh product brochure
- Phasix[®] Mesh IFU
- Characterization of host response, resorption, and strength properties, and performance in the presence of bacteria for fully absorbable biomaterials for soft tissue repair, Stoikes et al., 2017
- Poly-4-Hydroxybutyrate (P4HB) surgical Mesh for Hernia Repair, Gandhi et al.

Referenzen

- #1 Internal market research, data on file, 2015.
- #2 Liang MK, Li LT, Nguyen MT, Berger RL, Hicks SC, Kao LS. "Abdominal reoperation and mesh explantation following open ventral hernia repair with mesh." *The American Journal of Surgery* 208.4 (2014): 670-76.
- #3 Itani KM, et al. "Prospective study of single-stage repair of contaminated hernias using a biologic porcine tissue matrix: the RICH Study." *Surgery* 2012; 152(3): 498-505.
- #4 Annor AH, Tang ME, Pui CL, Ebersole GC, Frisella MM, Matthews BD, Deeken CR. "Effect of enzymatic degradation on the mechanical properties of biological scaffold materials." *Surgical Endoscopy* 26.10 (2012): 2767-778.
- #5 Harth KC, et al. "Effect of surgical wound classification on biologic graft performance in complex hernia repair: an experimental study." *Surgery* 2013; 153(4): 481-492.
- #6 Long-term clinical data not yet available.
- #7 Preclinical data on file at C. R. Bard, Inc. Results may not correlate to clinical performance in humans.
- #8 Nguyen PT, Asarias JR, Pierce LM. "Influence of a new monofilament polyester mesh on inflammation and matrix remodeling." *J Invest Surg* 2012; 25: 330
- #9 Bryan N, Ahswin H, Smart NJ, Bayon Y, Hunt JA. "In vitro activation of human leukocytes in response to contact with synthetic hernia meshes." *Clin Biochem* 2012; 45: 672
- #10 Aydinuraz K, Agalar C, Agalar F, Ceken S, Buruyurek N, Voral T. "In vitro S. epidermidis and S. aureus adherence to composite and lightweight polypropylene grafts." *J Surg Res* 2009; 157: e79.
- #11 Amid PK, Shulman AG, Lichtenstein IL, Hakaha M. "Biomaterials for abdominal wall hernia surgery and principles of their applications." *Langenbecks Archive Chir* 1994; 379(3): 168-71.
- #12 Klinge U, Junge B, Spellerberg B, Piroth C, Klosterhalfen B, Schumpelick V. "Do multifilament alloplastic meshes increase the infection rate? Analysis of the polymeric surface, the bacterial adherence, and the in vivo consequences in a rat model." *J Biomed Mater Res (Appl Biomater)* 2002; (63): 765-771.
- #13 Preclinical data on file at C. R. Bard, Inc. Phasix[®] ST Mesh shows similar mechanical strength over time to Phasix[®] Mesh in a preclinical model. Results may not correlate to clinical performance in humans.

- #14 Deeken CR, Matthews BD. "Characterization of the mechanical strength, resorption properties, and histologic characteristics of a fully absorbable material (Poly-4-hydroxybutyrate—Phasix[®] Mesh) in a porcine model of hernia repair." *ISRN Surgery* 2013; 2013: 238067. doi: 10.1155/2013/238067. PubMed PMID: 23781348; PubMed Central PMCID: PMC3679684.
- #15 Deeken CR, Abdo MS, Frisella MM, Matthews BD. "Physicomechanical evaluation of absorbable and nonabsorbable barrier composite meshes for laparoscopic ventral hernia repair." *Surgical Endoscopy* 25.5 (2010): 1541-552.
- #16 Estimated from Standard Curve in manuscript (Martin, et al. JSR 2013).
- #17 *Biochemical Engineering Journal*. 2003; (16): 97-105.
- #18 *International Journal of Polymer Science*. 2012; 1-11.
- #19 *Langenbecks Archives of Surgery*. 2012; (397): 363-371.
- #20 Blatnik A, Krpata D, Jacobs M, Gao Y, Novitsky Y, Rosen M. In Vivo Analysis of Morphologic Characteristics of Synthetic Mesh to Resist MRSA Adherence. *J Gastroint Surg*. 2012; 16(11): 2139-44.
- #21 Ceydeli A, Rucinski J, Wise L. "Finding the best abdominal closure: An evidence-based review of the literature." *Current Surg*. 2006; 62: 220-225.
- #22 Chu CC. The effect of pH on the in vitro degradation of poly(glycolide lactide) copolymer absorbable sutures. *J Biomed Mater Res*. 1982 Mar; 16(2): 117-24.
- #23 Katz et al. New synthetic monofilament absorbable suture made from polytrimethylene carbonate. *Surgery, Gynecology & Obstetrics*. 1985; 161(3): 213-222.
- #24 Novitsky YW. "Bridging versus closing the defect during laparoscopic ventral hernia repair." *The SAGES Manual of Hernia Repair* (2012): 439-44.
- #25 Gonzalez AM, et al. "Laparoscopic ventral hernia repair with primary closure versus no primary closure of the defect: Potential benefits of the robotic technology." *Int J Med Robotics Comput Assist Surg* 11.2 (2014): 120-25.
- #26 Bacterial clearance properties of a novel fully-absorbable barrier mesh (Phasix[®] ST) compared to a non-crosslinked acellular porcine dermal matrix (Strattice[®]) in a dorsal rabbit model. Corey R. Deeken, PhD*¹ & Jeffrey R. Scott, PhD*²; *¹Covalent Bio, LLC, St. Louis, MO; *²C. R. Bard, Inc. (Davol), Warwick, RI; DAV/PHST/0915/0019a

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Wenn Mängel am Produkt auftreten, wenden Sie sich wegen der Rücksendung unter 1-800-556-6275 an Davol, Inc. Copyright © 2010, 2013, 2014 C. R. Bard, Inc. Alle Rechte vorbehalten.
Bard, Davol, and Phasix sind Marken und/oder registrierte Marken von C. R. Bard, Inc. Copyright © 2016 C. R. Bard, Inc. Alle Rechte vorbehalten.
C. R. Bard GmbH · Wachhausstraße 6 · 76227 Karlsruhe
Bard Medica S.A. · Rinnböckstraße 3 · 1030 Wien
Bard Medica S.A. · Seestrasse 64 · 8942 Oberrieden/Zürich
Septra und Sepramesh sind Marken und/oder eingetragene Marken der Genzyme Corporation lizenziert an C. R. Bard, Inc.